


GUARDIAN SUCTION LINER

NEW ALPS™ Guardian Suction Liner - Features raised GripGel™ bands that grip the socket wall to form a secure interface between the socket and liner to prevent slippage or premature release. Low modulus GripGel™ bands stretch against the socket wall, while the inner wall conforms easily to the residual limb, to ensure there is NO RESTRICTION of blood flow or stiffening to inhibit donning.



Transtibial (BK) rings begin 3 1/2 inches from the center of the distal end.

Transfemoral (AK) rings begin 5 1/2 inches from the center of the distal end.

Key Features:

- Can be fully inverted for ease of donning
- Low modulus grip gel bands help ensure there is no restriction of blood flow
- No seams - single piece construction improves durability

Circle # 167 on Readers' Service Card



Use the Skin Reliever (ENCP) with the Guardian Liner to compensate for volume fluctuations with continued use of same socket.