



# GUARDIAN SUCTION LINER



- Full length inversion for ease of donning
- GripGel™ bands provide increased skin comfort
- High Performance Knitted Fabric Cover
- Transtibial (BK) rings begin 3 1/2 inches from the center of the distal end.
- Transfemoral (AK) rings begin 5 1/2 inches from the center of the distal end.

Circle # 167 on Readers' Service Card



In Partnership with:  
  
More than distribution.  
Partnership.

CE  
Tel: 727.528.8566 Tel: 800.574.5426  
www.easyliner.com info@easyliner.com

